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GB TABLE OF CONTENTS

1	For your safety and that of your patients	20
1.1	Symbols and information on the appliance.....	20
1.2	General safety information	20
1.3	Safety information for working Rules and guidelines ..	22
2	Product description	22
2.1	Intended purpose	22
2.2	Set-up and mode of operation ..	23
3	About the instructions for use.....	23
3.1	Scope.....	23
3.2	Target audience	23
3.3	Improvements	124
4	Transport and packaging	24
5	Operation	24
5.1	Commissioning	25
5.2	Decommissioning and storage... ..	26
6	Maintenance and care ..	27
6.1	Cleaning and disinfecting	28
6.2	Performance test	28
6.3	Maintenance	29
6.4	Warranty, customer service and dispatch	29
6.5	Troubleshooting	30
6.6	Disposal	30
7	Combination with downstream appliances and accessories	30
8	Appendix	31
8.1	Technical specifications.....	31
8.2	Flow curves (air)	32
8.3	Classification / UMDNS code.....	33

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1 FOR YOUR SAFETY AND THAT OF YOUR PATIENTS

1.1 Symbols and information on the appliance

GB

The following section offers an explanation of the meaning of the symbols that are both attached to the appliance (e.g. type plate) and contained in or on the accompanying information and packaging.

	<p>"Caution"</p> <p>This symbol indicates instructions which, in the event of non-compliance, could result in injury and/or property damage.</p>
	<p>"Beware"</p> <p>This symbol indicates instructions which, in the event of non-compliance, could result in property damage.</p>
	<p>"Free from oil and grease"</p> <p>This symbol denotes equipment that uses medical gases and must therefore be kept free from oil and grease.</p>
	<p>"Manufacturer"</p> <p>This symbol denotes the name and, possibly, address of the manufacturer.</p>
	<p>"Date of manufacture"</p> <p>This symbol denotes the date of manufacture.</p>
	<p>"Temperature restriction"</p> <p>This symbol indicates a temperature range, e.g. in connection with storage or operating conditions.</p>

1.2 General safety information

Please read and observe the following information before using or performing maintenance on the TESCOM EUROPE ZV compressed gas regulator II.

- The pressure regulator may only be used for its intended purpose as described in these instructions for use (see Section "Intended purpose").
- These instructions for use are integral to the appliance and must be available to the user at all times.
- You must understand how the appliance works before commencing operation with it.
- Never use an appliance if it is leaking or faulty!
- Medical gases are a medicinal product and may only be administered under the supervision of qualified medical staff who are familiar with the precautions that need to be observed for the respective treatment.
- Before commissioning the appliance, make absolutely sure that the gas type, the gas volume flow range and the connection are correct and that any upstream or downstream equipment has been fitted securely.
- Not suitable for operation in the range of strong magnetic fields, e.g. those generated by magnetic resonance imaging.
- Only accessories bearing the CE mark that are compatible with the appliance both in terms of intended use and connection may be used (see Section "Combination with downstream equipment").
- Never try to open the pressure regulator since this could result in the release of spring-loaded parts which might cause injury.
- To avoid the risk of infection or bacterial contamination, please read and observe the information in the Section "Maintenance and care".
- The pressure regulator must not be immersed in liquids or sterilised.
- Foreign matter must not get inside the appliance. This applies to dust and disinfectants in particular. Use the supplied protective caps / plugs.
- You **MUST** wash your hands before working on the oxygen supply.
- Ensure that the pressure regulator cannot come into contact with oil or grease during operation or storage. Oil and grease, as found for example in hand creams, can react violently with pressurised gases (especially oxygen and nitrous oxide). Risk of explosion!
- Do not use any lubricants.
- Smoking and naked flames are strictly forbidden in the vicinity of valves and fittings carrying oxygen.
- The connections on the appliance must match those on the gas supply or downstream equipment. Adapters must not be used.
- Only suitable for stationary operation in hospitals and clinics.
- Do not perform any technical or design modifications.
- Comply with the storage and operation conditions.
- Do not use the appliance for direct supply to the patient.
- Do not switch several appliances in sequence.

- Do not apply any external loads, support connected appliances.

- Do not hang objects on it.

- Protect against impacts, jolts, falls and vibrations.

1.3 Safety information for working with medical gases

GB

Oxygen (O₂) increases combustibility and flammability as well as encouraging flammable materials to burn. Greases and oils (such as those found in hand creams or lubricants) can ignite on contact with pressurised oxygen. Good ventilation is absolutely essential. High levels of concentration can be damaging to health. Smoking and naked flames are forbidden.

Nitrous oxide (N₂O) has a strong anaesthetic effect. It increases combustibility and flammability as well as encouraging flammable materials to burn. Greases and oils (such as those found in hand creams or lubricants) can ignite on contact with pressurised nitrous oxide. High levels of concentration could cause anoxia and apnoea. Good ventilation is absolutely essential. Smoking and naked flames are forbidden.

Carbon dioxide (CO₂) has an intoxicating effect at low concentrations, and a debilitating effect (apnoea) at higher concentrations. Good ventilation is absolutely essential.

High levels of **nitrogen (N₂)** can reduce O₂ content in breathable air in enclosed spaces. Ventilation is absolutely essential.

1.4 Rules and guidelines

Compliance with national laws, regulations and guidelines on the safe handling and operation of medical equipment in the context of 2007/47/EC and valves and fittings that carry gas, on handling gas under pressure, especially oxygen (O₂) and nitrous oxide (N₂O), with the information provided by the suppliers and manufacturers of the pressurised gases and with the instructions for using any connected equipment is mandatory.

2 PRODUCT DESCRIPTION

2.1 Intended purpose

The TESCOM EUROPE ZV compressed gas regulator II is used for a variety of medical equipment in the field of medicine and therapy to regulate the pressure at which medical gases are delivered.

- For direct connection by a trained user to an extraction point of a central gas supply system.
- Only for supplying suitable medical equipment in accordance with the technical specifications indicated on the type plate and in the instructions for use.
- Only for stationary use in hospitals and clinics.
- Only to be used for the medical gas indicated on the appliance.

Appliances with a NIST inlet connection may only be used with a TESCOM EUROPE central supply hose of up to 3 m in length that is connected to the inlet and used in configuration with a quick connector and NIST spigot nut.

The pressure regulator is not a shutoff valve.

2.2 Set-up and mode of operation

The TESCOM EUROPE ZV compressed gas regulator II is a single-stage, direct-acting, spring-loaded pressure regulator without primary pressure compensation. A metal membrane acts as the seal to the atmosphere. The downstream pressure is either fixed or can be adjusted by handwheel, depending on the model.

The version with universal connector (Air/O₂) at the inlet is fitted with an Air connector at the outlet.

3 ABOUT THE INSTRUCTIONS FOR USE

3.1 Scope

These instructions for use are valid for the TESCOM EUROPE ZV compressed gas regulator II with the part number indicated on the last page of these instructions for use.

Additional instructions for use must be observed for certain versions.

3.2 Target audience

- These instructions for use are intended for anyone responsible for the installation, operation or monitoring of the above valves and fittings.
 - The operator has a duty to ensure that only personnel who are familiar with the general practice and safety measures associated with working with gases and gas valves and fittings, as well as with the fundamental regulations with regard to occupational safety and accident prevention, are permitted to work on the device. These personnel must also have been taught how to handle the pressurised device, and have read and understood the instructions for use.
- The design of the appliance, the accompanying documentation and the marking are based on the assumption that the appliance will only be installed, operated and monitored by trained and qualified personnel and that the trained personnel are familiar with certain general characteristics of the appliance. As such, the instructions, warnings and safety information focus mainly on features that are particular to the TESCOM EUROPE appliance. The accompanying documentation does not address various general hazards and risks that should be obvious to qualified medical personnel. These include the incorrect use of the equipment and the possible side effects on patients in an unusual state of health.

- The users of appliances for medical gases are themselves responsible for taking the necessary safety precautions to ensure the safe operation of the appliances. A number of precautions can be taken to ensure patient safety. The user is responsible for deciding which precautions to take.
- It is absolutely essential that you read and understand at least the Sections "For your safety and that of your patients" and "Installation and operation" before attempting to operate the valves and fittings.
- Operators must be given access to these instructions. Additional copies can be supplied on request.
- Maintenance may only be performed by TESCOM EUROPE or companies authorised by TESCOM EUROPE.

3.3 Improvements


If you have any suggestions for improving the product or the instructions for use, we would be delighted to hear from you. Both our products and the instructions for use are subject to a process of continuous development and improvement. The address and telephone number for TESCOM EUROPE GMBH & CO. KG appear on the last page of these instructions for use.

4 TRANSPORT AND PACKAGING

The appliance may only be transported (e.g. delivery to the manufacturer) and stored in suitable, sturdy packaging (e.g. the original packaging). During transport and storage the inlet and outlet connections must be covered by suitable and clean protective caps/plugs.

5 OPERATION

ZV Compressed gas regulator II

Legend

- 1 Handwheel
- 2 Outlet pressure gauge
- 3 Inlet connection
- 4 Outlet connection
- 5 Type plate

- The user and/or operator are responsible for using and operating the appliance and accessories.
- You must wash your hands thoroughly before commencing any work involving the oxygen supply. Hydrocarbon compounds (e.g. oils, greases, rubbing alcohols, hand creams or remains of adhesive plasters) can cause explosive reactions if they come into contact with highly compressed gases.
- Medical gases are a medicinal product and as such may only be used in accordance with instructions issued by the doctor responsible for the care of the patient in question.

5.1 Commissioning

- If the appliance has been out of service for more than 6 months, a performance test must be performed before using it.

- All of the connections to the gas extraction point, the appliance, connected equipment and accessories must be undamaged, clean and free from oil and grease! Soiling on the pressure regulator can ignite.

- The gauge shows the outlet pressure relative to the ambient pressure. Observe the ambient pressure (large height, overpressure) when setting the pressure.

5.1.1 Adjustable pressure regulator

1. Check on the type plate whether the appliance is suitable for the intended purpose (type of gas, pressure). The maximum permissible inlet pressure of the appliance must be greater than or equal to the line pressure.

2. Turn the handwheel as far as it will go to the left (anti-clockwise).

Only use moderate force to turn the handwheel as far as the stop.

3. Prepare the downstream equipment and the accessories required for operation in accordance with the relevant instructions for use.
4. Remove the protective caps covering the connections on the appliance and keep them in a safe place.
5. Appliances with NIST inlet connection: Connect the TESCOM EUROPE hose with inlet connection in accordance with the relevant instructions for use.

The temperature of the gas and the pressure regulator drops during operation. This can produce hypothermia in the patient.

6. Connect the inlet connection on the pressure regulator to the gas extraction point.
7. Use the correct and proper procedures to connect the downstream equipment and accessories to the outlet connection on the appliance.

The connections must match. Do not use adapters nor modify the connections on the appliances.

8. Read the outlet pressure on the gauge and, if necessary, adjust by turning the handwheel.

9. Do not leave the pressure regulator and the patient unsupervised during operation.

GB

To increase the outlet pressure

To reduce the outlet pressure (only while gas is flowing)

5.1.2 Non-adjustable pressure regulators

1. Check on the type plate whether the appliance is suitable for the intended purpose (type of gas, pressure). The maximum permissible inlet pressure of the appliance must be greater than or equal to the line pressure.

2. Prepare the downstream equipment and accessories in accordance with the relevant instructions for use and use the correct and proper procedures to connect them to the outlet connection on the appliance. The connections must match. Do not use adapters nor modify the connections on the appliances.

3. Appliances with NIST inlet connection: Connect the TESCOM EUROPE hose with inlet connection in accordance with the relevant instructions for use.

4. Remove the protective caps covering the connections on the appliance and keep them in a safe place.

5. Connect the inlet connection on the pressure regulator to the gas extraction point.

The temperature of the gas and the pressure regulator drops during operation. This can produce hypothermia in the patient.

Observe the minimum operating temperature of the downstream equipment and accessories. The temperature of the gas and the pressure regulator drops during operation.

5.2 Decommissioning and storage

5.2.1 Adjustable pressure regulators

Turn the handwheel on the pressure regulator as far as it will go to the left (anti-clockwise) and allow the gas in the pressure regulator to escape.

Only use moderate force to turn the handwheel as far as the stop.

1. Disconnect the downstream equipment and accessories from the appliance.
2. Disconnect the appliance from the gas extraction point.
3. Fasten the protective caps on the connections.
4. Store the pressure regulator in a clean, dust-free and dry place until the next time it needs to be used.

5. If the pressure regulator is out of use for a prolonged period, it should be cleaned (see Section "Maintenance and care").

It is also vital to ensure that maintenance intervals are observed for appliances in storage. Only pressure regulators that have been properly maintained whilst in storage may be put to use when taken out of storage.

5.2.2 Non-adjustable pressure regulators

1. Disconnect the pressure regulator from the gas extraction point.
2. Disconnect the accessories from the appliance.
3. Fasten the protective caps on the appliance's connections.

4. Store the pressure regulator in a clean, dust-free and dry place until the next time it needs to be used.

5. If the pressure regulator is out of use for a prolonged period, it should be cleaned (see Section "Maintenance and care").

It is also vital to ensure that maintenance intervals are observed for appliances in storage. Only pressure regulators that have been properly maintained whilst in storage may be put to use when taken out of storage.

6 MAINTENANCE AND CARE

Non-compliance with the maintenance and care instructions for the appliance can result in severe injury to the patient.

Foreign matter must not get inside the appliance. The pressure regulator may not be immersed in liquid!

The pressure regulator may not be sterilised!

6.1 Cleaning and disinfecting

- All of the inlet and outlet connections on the appliance must be sealed with suitable and clean protective caps/plugs prior to cleaning and disinfecting the equipment.

6.1.1 Cleaning

Wipe the surface of the pressure regulator with a dry disposable cloth or with a disposable cloth that has been moistened with clean water.

GB

Cleaning agents containing alcohol or emollients may ignite if they come into contact with compressed oxygen or nitrous oxide.

6.1.2 Disinfection

Only use approved surface disinfectants. Observe the manufacturer's instructions for use.

Disinfectants that are based on the following materials are compatible:

- aldehydes
- quaternary ammonium compounds.

Disinfectants based on the following materials are not suitable as they may damage the materials:

- halogen-splitting compounds
- strongly organic acids
- oxygen-splitting compounds

Users in the Federal Republic of Germany are recommended to always use disinfectants that are included in the latest DGHM register (DGHM: Deutsche Gesellschaft für Hygiene und Mikrobiologie, mhp-Verlag GmbH Wiesbaden). The DGHM register also lists the active agents on which each disinfectant is based.

The recommendation based on the active agents listed above applies in those countries that are unaware of the DGHM list.

6.2 Performance test

A visual inspection must be performed at least every 6 months in brightly lit, glare-free conditions to check for mechanical damages. The following must be inspected:

- Inlet connection
- Housing
- Pressure gauge
- Handwheel
- Outlet connection
- Legibility of type plate and gas type plate
- Alignment of connections and gauge to the appliance axes
- Secure fit of the connections and gauge

The handwheel on adjustable appliances must turn easily.

It is also vital to ensure that maintenance intervals are observed for appliances in storage. Only pressure regulators that have been properly maintained whilst in storage may be put to use when taken out of storage.

Should faults be detected during performance testing or operation, the pressure regulator must be taken out of service immediately. Contact TESCOM EUROPE or an authorised specialist dealer to have the fault rectified.

6.3 Maintenance

Maintenance including - amongst other things - the replacement of all wear parts must be carried out every 6 years. Non-compliance with this deadline will invalidate the warranty provided by TESCOM EUROPE.

The safety and function of the appliance can only be assured subject to maintenance being carried out by TESCOM EUROPE or an authorised company and provided that original TESCOM EUROPE spare parts are used.

6.4 Warranty, customer service and dispatch

All products purchased from TESCOM EUROPE GMBH & CO. KG are covered by the statutory warranty from the date of purchase.

The warranty obligations and liability of TESCOM EUROPE GMBH & CO. KG for any consequences are excluded in the following cases:

- Non-compliance with the instructions for use
- Use for any other than the intended purpose
- Operating errors
- Improper use or handling
- Negligence
- Tampering with the pressure regulator by unauthorised third parties, e.g. opening the appliance or replacing connection pieces
- Act of God
- Damage during transport caused by improper packaging
- Failure to maintain the appliance
- Failure to use original spare parts

The warranty and liability conditions of the conditions of sale and delivery of TESCOM EUROPE GMBH & CO. KG are not extended by the above notices.

Please contact our repairs department before returning any goods for repair or any faulty goods. In the interests of our staff we also require a declaration indicating whether – and if so which – substances have been used with the pressure regulator and which measures (e.g. flushing) have been taken so that we can exclude any hazards during repair.

You will find the manufacturer's address on the last page of these instructions for use.

6.5 Troubleshooting

GB

Problem	Cause	Solution
Connection between gas extraction point and appliance is not leak-tight.	Inlet connection is damaged. Gasket on the gas extraction point is damaged.	Must be repaired by TESCOM EUROPE or an authorised specialist dealer.
Outlet pressure increases	Soiled or damaged valve seat.	Must be repaired by TESCOM EUROPE or an authorised specialist dealer.
Leakage in the housing area	Faulty diaphragm.	Must be repaired by TESCOM EUROPE or an authorised specialist dealer.

6.6 Disposal

Please contact an approved and certified disposal company to ensure the proper disposal of the appliance. Your environmental protection representative or municipal authorities will be able to provide the relevant contact details.

7 COMBINATION WITH DOWNSTREAM APPLIANCES AND ACCESSORIES

The user and/or operator are responsible for using and operating the pressure regulator and the equipment connected to it.

Downstream appliances and accessories must comply with the following criteria:

- It must be a medical product with CE mark.
- It must be marked with the same type of gas as the pressure regulator.

Risk of confusing gas types when using non-gas-type-specific connections (e.g. hose nozzle).

Risk of the hose slipping if hose nozzles are used.

- Downstream appliances and accessories must be suitable for the pressure regulator and the application. Details about gas type, outlet pressure and flow can be found in the technical specifications in the appendix to these instructions for use.

The temperature of the gas and the pressure regulator drops during operation. This can produce hypothermia in the patient.

Should you have any questions or require specialist advice, please contact us or visit us on the Internet at www.tescom-europe.com.

You will find the manufacturer's address on the last page of these instructions for use.

8 APPENDIX

8.1 Technical specifications

Type of gas:O₂, Air, CO₂, N₂O

Permissible inlet pressure p_{max}:

D45715 / D513704 + 1 bar

D45714 / D513698 + 2 bar

Maximum outlet pressure p₂:

D45715 / D513704 bar

D45714 / D513698 bar

Nominal gas flow Q₁:

D45715 / D5137020 l/min

D45714 / D5136925 l/min

The actual flow is dependent on the pressure difference in the pressure regulator (see Section "Flow curves").

Permissible ambient conditions

Air humidity:5% to 95% rel. humidity,
non-condensing



Air pressure

during operation:900 hPa to 1100 hPa

Storage:115 hPa to 2000 hPa

Temperature
Storage


Operation


Weight:650 g

8.2 Flow curves (air)

D45715 / D51370

8.3 Classification / UMDNS code

Pursuant to EC Directive 2007/47/EC

Appendix IX Class: IIb

Nomenclature for medical devices 13-320 compressed gas regulators

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